

Network Security with Mettler-Toledo

EasyMax[®]/OptiMax[™] Instruments

This document provides an overview of measures implemented on Mettler-Toledo EasyMax[®]/OptiMax[™] instruments to guard against security threats resulting from malware or network attacks.

Overview of EasyMax[®]/OptiMax[™] hardware and software:

Item	Description
Operating System	Microsoft Windows CE 5.0
Instrument Control Software	Custom Mettler-Toledo SW components running on .NET Compact Framework 3.5 and C/C++.
Processor Board	Toradex Colibri with Marvell PXA320 Processor
Ports for external connections	<ul style="list-style-type: none"> • USB Used to connect USB sticks to export experiment data or to upgrade the instrument firmware. • Ethernet / RJ45 Used to connect the instrument to the local network. PC software may be used to control or service the instrument.

Mettler-Toledo AG, AutoChem RXE
 Sonnenbergstrasse 74
 CH-8603 Schwerzenbach, Switzerland
 Phone +41-44 806 77 11
 Fax +41-44 806 72 90

Internet www.mt.com/autochem
 E-Mail icare@mt.com

© 03/2012 Mettler-Toledo AG

Malware on Windows[®] CE

Unlike Windows[®] Desktop Operating Systems (OS) such as XP, Vista or Windows[®] 7, malware is extremely rare on Windows CE due to the following reasons:

- Most malware is specifically developed for Desktop OS and does not run on Windows[®] CE.
- The tools that are typically used as entry points for malware, such as E-Mail readers, PDF readers, and Internet browsers, are not available on Mettler-Toledo instruments.
- There is no built-in OS auto-start mechanism for executables in the top folder of USB sticks.

Protection against network attacks

The following measures were implemented in order to protect the instruments from network attacks:

- All incoming ports except the ones needed for communication with PC software (TCP 8001, 8018, and 8019) are locked down (firewalled).
- In the event an instrument is flooded with IP packages (Denial of Service Attack), the service that answers network requests will be shut down and the instrument will transition to a safe temperature automatically.

Safe handling of user credentials

No user authentication is required to work with EasyMax[®]/OptiMax[™] instruments. Therefore, no credentials, such as passwords, are stored on the instrument or transmitted over Ethernet.