

# iC PVM™ Quick Start Guide

This guide introduces a new user to the iC PVM graphical user interface (GUI) and describes several choices for getting started. Create and record an experiment to save images that best represent a particle system and its process dynamics. Each experiment can easily generate a report, image slideshow, or video.

## Start Page

A start page appears when you first open iC PVM. At a glance you can see instrument status and begin an experiment, open a recent experiment file, or import an image sequence file (\*.seq).

**Click/tap Quick Start to quickly create a new iC PVM experiment.**

**Open recent experiments.**

**Browse to other experiment files**

**Import a sequence file**

**Open links to:**

- Activate software licenses
- View instrument configuration
- Check options including templates and iC Data Center
- Access iC PVM information and user assistance
- Go to AutoChem customer website

**Recent**

- Exp 2015-01-13 10-43  
Images: 39  
CVC PVM image files
- Exp 2015-01-09 10-24  
Images: 49  
CVC PVM image files
- Exp 2015-01-12 10-00  
Images: 1  
CVC PVM image files

Open other files ...

**Links**

- Licensing
- Instrument
- Options
- ? Help
- Customer Community Site

**New**

- Quick start**  
Create an experiment with default settings
- Clone Experiment**  
Use settings from a previous experiment
- From Template**  
Use settings from a user defined template

**Video Learning**

- Take a Tour**  
Introduction to iC PVM  
See a brief overview of the software
- Ensuring Quality Images**  
Learn about focusing and lighting
- Sharing Data and Images**  
Report your results in useful formats

**Glance at instrument connection status.**

**Images appear when instrument is connected, powered, and inserted in a particle/droplet system. NOTE: A clean probe window appears black when it is not inserted in a particle system.**

**Links to product video tutorials.**

**Alternatives for starting a new experiment:**

- Cloning an existing experiment (tweak settings)
- From a template (advanced option)

**Next**—Name experiment file and set timing options...

## Name Experiment File and Set Basic Options

After you click/tap **Quick Start** on the Start Page, click **Create** to accept the default filename and location of the experiment file that will have either the default or specified save interval and duration settings.

Edit or accept experiment file name and location.  
If using iC Data Center, name the experiment according to your Standard Operating Procedure.

You can change these interval and duration settings now or during the experiment.

The 'Quick Start' dialog box contains the following sections:

- Create an experiment with default settings**: Indicated by a sun icon.
- File Specifications**:
  - Experiment Name: Exp 2015-01-07 10-01
  - Folder: C:\iC PVM image files (with a browse button '...')
- Experiment Duration**:
  - Save Interval: 10 seconds (with a dropdown arrow)
  - Experiment Duration: 8 hours (with a dropdown arrow)
- Create**: A large blue button with a sun icon and the text 'Create'.

Amount of time between each image saved.

Experiment file holds up to 6,000 images.

Click/tap **Create** to make the experiment.

**Next**—Begin  
Recording images to  
the experiment file...




## Record Images

After you click/tap **Create** to make the experiment, the main iC PVM workspace appears with the LIVE ribbon. NOTE: Live images on display are unsaved until you click/tap the record button (●). When recording starts, the LIVE experiment toolbar changes to green.

Adjust, as needed: Focus (V19 only), Lighting, Save interval (and experiment phases)

Click/tap (●) to start recording images.




Multiple **ribbon tabs** show options and tools for the iC PVM workspace. (LIVE ribbon for ParticleView V19 shown)

LIVE experiment toolbar

Begin 'pinning' selected images by:

- Clicking/tapping the Pin button (from either location), or
- Entering a process note



HOME ribbon shown

Partial Zoom

**Next**—The main workspace shows images only, by default. But a split view layout includes Images and Trends...




## Observe Images and Trends

iC PVM includes a trends viewer and filmstrip of the saved images in sequence. The sample below shows a split layout of the two main sections in the iC PVM workspace—**Images** at the top and **Filmstrip and Trends** on the bottom.

To change the layout:

- Click/tap the split icon ( ). Notice the icon changes to the full version ( ) on the Images bar.
- Or, select **Image and Trends** layout from the VIEW ribbon.




Import an iControl or iC FBRM trend.

HOME ribbon shown

This is the **Image and Trends** layout where the workspace layout is split into two main sections:

1. **Images** collected

2. **Filmstrip** of saved images, and **Trends** chart of image trends such as Relative Backscatter Index (RBI) or trends imported from iC FBRM or iControl



Notes panel

- When you pin an image, it appears here.
- If panel is collapsed, click/tap Notes bar to expand it.

Enter note on pinned images.

Double-click/tap in an image thumbnail to move back and forth between the main workspace and the **Pinned images** view, described next...

**Next**—Review all the images pinned from the experiment...





## View Pinned Images

When you click/tap the **Pinned** button (  ) in the Notes panel or from the VIEW ribbon, the iC PVM workspace displays a gallery of all pinned images in this experiment.

**VIEW ribbon shown**

Select preferred image size

**To return to main workspace:**

- Double-click/tap on an image, or
- Click/tap the view button, or
- Click/tap one of the first two layout buttons in the VIEW ribbon.


Expand/collapse image properties including notes, image analysis, trend data, and acquisition information for a selected image.

Enter/edit notes (in either place).

Review image analysis data.

Expand/collapse technical details about image acquisition.

For completed experiments, iC PVM includes automatic reporting and a choice of export options (slide show, video, sequence files). Access both options from the HOME ribbon and FILE menu.



**VIEW ribbon:**

- FILE
- HOME
- IMAGE
- VIEW (Active)
- LIVE

**VIEW ribbon options:**

- Image Only
- Image and Trends
- Pinned Images
- Add Viewer
- Rulers
- Mini Scale
- Small images
- Medium images
- Large images
- Compress data gaps
- Show trends
- Show filmstrip
- Show all axes
- Pinned images only
- Filmstrip and Trend Viewer
- Time Format
- Local
- Elapsed

**Image details panel:**

**Selected image:** 4 of 4 pinned images

**Experiment name:** Exp 2015-01-07 10-01

**Local time:** 01/07/2015 1:46:38.66 PM

**Elapsed time:** 03:45:17.16


**Note:** Endpoint-Fine crystal formation

**Image analysis:**

- Total effective laser power: 0.00
- Relative backscatter index: 140.88
- Information content: 27.18 %
- Avg image brightness: 55.25 %
- Gray level @ 96.0 pcntile: 67.07 % gray

**Acquisition info:**

- Pixel width: 1360
- Pixel height: 1024
- Binned: False
- Accepted by: Accept all
- Laser state: 1 1 1 0 1 1 1 0
- Laser type: F F F B F F F B
- Laser intensity: 0.00
- Camera gain: 0.00
- Pixels/micron X: 1.25
- Pixels/micron Y: 1.25
- Focus pos. (um): 0
- Date Time: 1/7/2015 13:46:38.666
- Elapsed Time: 03:45:17.156



**IMAGE ribbon:** Apply common image adjustments to enhance an image, and mark up images using drawing/measuring tools.

**ANALYZE ribbon:** In completed experiments, you can set image acceptance properties as an aid in filtering saved images.



### Mettler-Toledo AutoChem, Inc.

7075 Samuel Morse Drive  
Columbia, MD 21046 USA  
Telephone +1 410 910 8500  
Subject to technical changes

© 1/2015 MK-PB-0118-AC Rev A DCN 2527

**METTLER TOLEDO**

[www.mt.com/pvm](http://www.mt.com/pvm)